



Acidex VGB Iodine 30

Antiseptic Solution 10% Povidone-Iodine Solution

PRODUCT PROPERTIES

As having a wide spectrum bactericidal effect, it's effective against gram (+) and gram (-) bacteria, fungi and viruses.

Especially used for treatment of contaminated wounds and antiseptics of skin mucous membrane.

RANGE OF APPLICATION

Recommended for antiseptic wound treatment, burns, mycotic and bacterial origin wounds, super infected dermatosis, hygienic and surgical disinfection.

APPLICATION

Use Directly.

For external use only without dilution.

Should not be used on people who are allergic to iodine.



ANTIMICROBIAL PROPERTIES

Staphylococcus aureus,
Candida albicans, Pseudomonas
aeruginosa, Escheria Coli

COMPOSITION

100 g Acidex VGB IODINE 30
contains;

10 g PVP-Iodine (1% w/w
available iodine)

Glycerine, ph buffer, deionised
water.

Product Code : GD 81-025
Weight/Volume : 250 ml
Packaging : 36 pcs.....

Animal Care

Antiseptic Wound Care Products



Fastec **Oxy Blue**

4,5 g Oxytetracycline

TOPICAL SPRAY

PRODUCT PROPERTIES

Oxytetracycline belongs to the group of tetracyclines and acts bacteriostatically against many Gram-positive and Gram-negative bacteria like Bordetella, Campylobacter, Chlamydia, E. coli, Haemophilus, Mycoplasma, Pasteurella, Rickettsia, Salmonella, Staphylococcus and Streptococcus spp. The action of oxytetracycline is based on inhibition of bacterial protein synthesis.

RANGE OF APPLICATION

An aid in the treatment of foot rot in sheep, digital dermatitis in cattle and superficial skin infections caused by oxytetracycline sensitive organisms in pigs, sheep and cattle.

APPLICATION

For topical use only. Shake well before use. Remove any dead tissue, dried serum or secretions and debris from infection site. Hold the spray can upright approximately 15-20 cm from the surface and cover the area completely once or twice daily and repeat for 3-5 days, or as required. The blue marker dye scours out completely from wool. Therefore, product is suitable for use on fibre-producing animals. Replace cap when not in use.



COMPOSITION

200 ml Fastec OXY BLUE

contains;

4.5 g Oxytetracycline Hydrochloride, Isopropyl Alcohol, Hydrocarbon gas, Dye, Deionized water.

Product Code : FV 007-02
Weight/Volume : 200 ml
Packaging : 36 pcs.