

# Personal Hygiene

## Products



## Fastec **Acihand**

### Antibacterial Foamy Hand Soap

#### PRODUCT PROPERTIES

- Contains antiseptic agents. Provides hygienic hand cleaning.
- Contains moisturizers, does not irritate the skin.
- Five times more economical than regular hand soap: while 1000 ml of regular liquid hand soap provides 300 washes, the same amount of Foamy Soap provides upto 1500 hand washes.

USAGE: Wet your hands with running water. Get 1-2 pumps of product on palm. Rub and lather. Scrub all

surfaces, including the backs of your hands, wrists, between your fingers and under your fingernails for 1-3 min. Rinse well. Dry hands with a sterile or disposable towel or air dryer. Use a sterile towel or your elbow to turn off the faucet.

#### ANTIMICROBIAL PROPERTIES

- Gram (+) and gram (-) bacteria, yeasts, molds, fungi (Staphylococcus aureus, Pseudomonas aeruginosa, Candida albicans, Escherichia coli).

#### COMPOSITION

Fastec ACIHAND contains;

- Sodium Lauryl Ether Sulphate
- Cocoamidopropyl Betaine & Glycerine
- 0,3 g Alkyl Dimethyl Benzyl Ammonium Chloride and deionised water.

Product Code	: FK 106-1	FK 106-5	FK 106-10	FK 106-20	FK 106-200
Weight/Volume	: 1 L	5 L	10 L	20 L	200 L
Packaging	: 16 pcs.	4 pcs.	2 pcs.	1 pcs.	1 pcs.

## Fastec **Acihand H 01**

### Hygienic & Surgical Hand Disinfectant/ Scrub

#### PRODUCT PROPERTIES

Fastec Acihand H 01 provides hand sanitation with its antibacterial effect. Due to moisturizers it contains and pH compatible with skin, doesn't irritate hands.

USAGE: Wet your hands with running water.

Pour 6-10 ml the product on palm directly or using the pump. Rub and lather. Scrub all surfaces, including the backs of your hands, wrists, between your fingers and under your fingernails

for hygienic hand wash: 2 x 3 ml - rub for 1 min.

for surgical hand wash: 2 x 5 ml - rub for 3 min.

Rinse well. Dry hands with a sterile or disposable

towel or air dryer. Use a sterile towel or your elbow to turn off the faucet.

#### ANTIMICROBIAL PROPERTIES

- Gram (+) and gram (-) bacteria, yeasts, molds, fungi (Staphylococcus aureus, Pseudomonas aeruginosa, Candida albicans, Escherichia coli).



#### COMPOSITION

Fastec ACIHAND H01 contains;

- 4 g Chlorhexidinedigluconat, glycerine, hydroxyethylcellulose and deionised water.



Product Code	: GD 28-1	GD 28-5	GD 28-10	GD 28-20	GD 28-200
Weight/Volume	: 1 L	5 L	10 L	20 L	200 L
Packaging	: 16 pcs.	4 pcs.	2 pcs.	1 pcs.	1 pcs.



## Acidex **Acisept Dz 21** Alcohol Based Hand Disinfectant

**RANGE OF APPLICATION:** It is an alcohol based hand disinfectant which provides effective hand hygiene in the medical settings such as hospitals, clinics, laboratories and ambulances. It can also be used as a surgical hand disinfectants.

**PRODUCT PROPERTIES**

- Bactericide, virucide, fungicide.
- Contains ingredients that protect and moisturize hands.
- Inactivates microorganisms within 30 seconds of application.

**USAGE:** Pour 3-5 ml of the product on palm directly or

using the pump. Let the product dry by rubbing into hands and scrubbing all surfaces, including the backs of your hands, wrists, between your fingers and under your fingernails for at least 30 seconds. Don't rinse.

**ANTIMICROBIAL PROPERTIES:** Effective against Hepatitis B, HCV / HIV, gram (+) and gram (-) bacteria (including TBC), fungi, SARS virus, inactivates Adeno, Polio and Rota viruses. (Staphylococcus Epidermidis ATCC 13047, Pseudomonas Aeruginosa ATCC 12228, Candida Albicans, Escherichia Coli ATCC 13047, Aspergillus Flavus)

Product Code	: FK107-1	FK107-5	FK107-10	FK107-20	FK107-200
Weight/Volume	: 1 L	5 L	10 L	20 L	200 L
Packaging	: 16 pcs.	4 pcs.	2 pcs.	1 pcs.	1 pcs.

**COMPOSITION**

Each 100 g of Solution Contains:  
70 g Isoprophyl Alcohol (CAS No: 67-63-0), Humectants, pH modifiers and deionized water.

## Acidex **Dr. Faber Hygiene** Alcohol Free Hand Disinfectant

**PRODUCT PROPERTIES:**

- Uniquely formulated with quaternary ammonium compounds as ingredients, which disintegrates bacterial cell membranes.
- Its updated, high tech renovated composition is effective against a wide range of bacteria and viruses.
- Economic due to foam form, provides 1500 washes.
- Gentle on the skin due to its skin regenerative, moisturizer and protective ingredients.
- Provides effective and quick disinfection.
- Phenol-free and environmentally friendly.

**RANGE OF APPLICATION:**  
Used in hospitals, clinics, laboratories, and food and beverage producing facilities.

**ANTIMICROBIAL PROPERTIES:** Effective against a wide range of bacteria and viruses, such as Noro, SARS, BVDV, Vaccinia, Rota, HBV, HIV, HCV, Adeno, fungal infections, Gram (+) and Gram (-) bacteria.



**COMPOSITION**

100 g of Acidex Dr. Faber Hygiene contains;  
• 0,3 g Alkyldimethylbenzylammonium chloride, glycerine and deionised water.

Product Code	: FK111-1	FK111-5	FK111-10	FK111-20	FK111-200
Weight/Volume	: 1 L	5 L	10 L	20 L	200 L
Packaging	: 16 pcs.	4 pcs.	2 pcs.	1 pcs.	1 pcs.

